

How To Raise A Nurse Call Alert

For use during a medical emergency where assistance is urgently required from on site emergency response personnel.



1) To raise an alert, press the **CALL** button in the middle

2) A **BLUE** light next to **SENT** and a single three second beep indicates your alert has been raised



3) The alert is pushed through to the TCU, smartphone app, and via SMS (subject to system design and software subscription)



4) Remain calm while you wait for personnel to respond

5) A **GREEN** light next to **RECEIVED** and three x one second beeps indicates someone is responding to your alert



6) If no one responds within three minutes, the **BLUE** light will go off and you will need to press the **CALL** button again

7) If you need to cancel the alert at any time, press the **CANCEL** button